


Percuplasty™

Percutaneous 3D Bunion Correction®




TREACE
Medical Concepts, Inc.
A Step Ahead.™

Percuplasty™

Percutaneous 3D Bunion Correction®

- **Reproducible Corrections**
Elegantly instrumented to reproducibly achieve 3D MIS corrections
- **Reliable Throws**
Large diameter cobalt chrome wires designed to reduce skiving
- **Rapid Insertion**
Self-drilling screws to allow for efficient placement



Ordering Information

3.0mm Screws		4.0mm Screws	
SF30XX (XX = Length) Available Lengths: 18mm – 50mm In 2mm increments		SF40XX (XX = Length) Available Lengths: 24mm – 64mm In 2mm increments	
Shannon Burr	Shannon Burr	Wedge Burr	
TM12013 2mm x 13mm	TM12220 2.2mm x 20mm	TM22913 3mm	

Before use of the system, the surgeon should refer to the appropriate instructions for use and surgical technique for complete warnings, precautions, indications, contraindications, and adverse events. Risks include, but are not limited to: infection, pain, discomfort from the presence of the implant, loosening of the implant, and loss of correction with nonunion or malunion. If any of these occur, additional treatments may be needed. Additional information about risks, warnings, and instructions is available at Lapiplasty.com/surgeons/labeling.



**To learn more,
visit Treace.com**